



 **SAMPLIX<sup>®</sup>**

Technology from the  
Heart of Europe

Evacuated Blood Collection Tubes

## Safe and simple blood collection with Samplix®

Samplix® is an evacuated blood collection system for the collection of human samples. It was developed by the product experts at Greiner Bio-One, one of the world's leading medical technology companies, which is based in Austria, right at the 'heart of Europe'.

The Samplix® tubes are available with a variety of additives, which prepare the samples for the particular analyses. This allows an exceptionally safe sampling, because a manual addition of the additive is no longer necessary.

### Quality starts with Samplix®

The precisely pre-dosed vacuum in our Samplix® Evacuated Blood Collection Tubes offers you the great and essential advantage of valid test results over conventional blood collection tubes:

The closed vacuum system significantly reduces a potential contamination with the sample material. All materials used in tube manufacturing are considered chemically resistant and inert. Due to the production of plastic (PET), Samplix® tubes offer many other advantages:

- + **minimal contact activation for minimal to no interference on assays**
- + **highly transparent for good visibility of sample integrity and quality**
- + **nearly unbreakable for low injury risk and sample loss**
- + **Low risk of contamination thanks to closed system**
- + **Compliant with the tolerance range (+/- 10 %) according to ISO 6710**
- + **Correct blood-additive mixing ratio ensures valid test results**

### Ordering information Samplix® Evacuated Blood Collection Tubes

Item n°	Description	Cap Colour	Vol.	Size	Rack	Carton
554101	Samplix® Tube CAT Serum Clot Activator	● Red	4ml	13x75	100 pcs.	1.600 pcs.
554201	Samplix® Tube CAT Serum Sep Clot Activator	● Dark Yellow	4ml	13x75	100 pcs.	1.600 pcs.
554301	Samplix® Tube LH Lithium Heparin	● Green	4ml	13x75	100 pcs.	1.600 pcs.
554401	Samplix® Tube LH Lithium Heparin Sep	● Light Green	4ml	13x75	100 pcs.	1.600 pcs.
554501	Samplix® Tube K2E K2EDTA	● Lavender	2ml	13x75	100 pcs.	1.600 pcs.
554502	Samplix® Tube K2E K2EDTA	● Lavender	3ml	13x75	100 pcs.	1.600 pcs.
554511	Samplix® Tube K3E K3EDTA	● Lavender	2ml	13x75	100 pcs.	1.600 pcs.
554512	Samplix® Tube K3E K3EDTA	● Lavender	3ml	13x75	100 pcs.	1.600 pcs.
556101	Samplix® Tube CAT Serum Clot Activator	● Red	6ml	13x100	100 pcs.	1.600 pcs.
556201	Samplix® Tube CAT Serum Sep Clot Activator	● Dark Yellow	5ml	13x100	100 pcs.	1.600 pcs.
556301	Samplix® Tube LH Lithium Heparin	● Green	6ml	13x100	100 pcs.	1.600 pcs.

The Samplix® Evacuated Blood Collection Tubes are used to collect, transport, store and process blood for testing serum, plasma or whole blood in clinical laboratory and are intended for professional and single-use.

**STORAGE**  
4-25°C (40-77°F)  
avoid exposure to direct sunlight

**SHELF-LIFE**  
12 months





Greiner Bio-One with headquarters in Kremsmünster, Austria, is a family company and a market leader in the fields of preanalytics, life sciences, diagnostics and OEM.

Our company stands for high quality. This can be seen in our innovative product solutions: based on our many years of experience as well as continual development, Greiner Bio-One ensures safety and reliability in daily routine work with patients in your healthcare facility.

We are committed to fulfilling your desire for health, joy and a long life for you and your family. Because Greiner Bio-One is YOUR **POWER FOR HEALTH!**

The Samplix® Blood Collection Tubes are manufactured at our state-of-the-art production plant in Thailand.

**FOR FURTHER INFORMATION  
PLEASE CONTACT US.**

**THAILAND**

Greiner Bio-One Thailand Ltd.  
**PHONE** +66 38 4656 33  
**FAX** +66 38 4656 36  
**E-MAIL** [office@th.gbo.com](mailto:office@th.gbo.com)